

Flexibility and comfort – an ideal solution for Acute and Residential care



The Invacare Jasmine Mobile lifter is a flexible solution that offers a high level of comfort for all clients. A versatile lifter, the Jasmine provides safe patient handling up to a 32 stone (200kg) weight capacity. With an expansive, wide-opening base and high lifting mechanism, the Jasmine offers easy manoeuverability and positioning for a wide range of lifts. With the latest control box technology, the Jasmine offers a 'Soft Start' and 'Soft Stop' function that provides a more comfortable and dignified transfer.

## Features and Options





#### **Spreader bar options** Scale

The **Jasmine** can be ordered with 2 or 4-point spreader bars. In widths: 450 mm and 550 mm.

A scale is available. Approved according to medical standard EN 45501, certified and sealed.



#### Linak<sup>®</sup> Jumbo<sup>™</sup> Care

The intelligent monitoring system includes a service reminder, emergency operation, audio battery alarm, and battery indicator.

#### Safety first, with an intelligent control box

The **Jasmine** is fitted with a control box that has an integrated system: an intelligent service lifetime monitoring system that highlights the battery status and service condition of the actuator to ensure the correct use and safety of the hoist.

### **Flexibility and optimized** maneuverability

The Invacare Jasmine is exceptionally mobile and easy to transport. The wide-opening base allows easy access for a wheelchair or armchair, whilst the high lifting mechanism is essential for transfers to and from elevated surfaces. Flexible and with optimised maneuverability, the Jasmine covers the needs of both the carer and client.

# **Technical data**

For more comprehensive pre-sales information about this product, including the product's user manual, please visit www.invacare.co.uk.



Jasmine - UK - 09/2017 - **1539003**