



rea

Confidence in posture solutions

rea

Azalea® Minor

Designed for petite adults
and teenagers

INVACARE
Yes, you can.

The **Rea Azalea Minor** is designed for those who require a seat that has a smaller depth and width. Adapted from the **Rea Azalea**, the **Azalea Minor** boasts all the advantages of a reliable, tilt-able wheelchair with its special weight-shifting mechanism. This model is ideal for all people who require a solid, stable wheelchair. The Azalea Minor is available with transit or self propel wheels, five frame colours and three upholstery choices.

Rea Azalea compatibility continues

Rea Azalea Minor designed for those who need a smaller-sized seat, such as petite adults or teenagers. Key features of **Azalea Minor** is the tilt-in-space weight-shifting mechanism, to ensure a stable tilt and makes manoeuvres easier. The width and depth of the **Azalea Minor** can be adjusted to suit individual needs.

Tension adjustable backrest

The tension adjustable backrest offers individual shaping to accommodate more complex spine shapes. This back style is available with a choice of positioning backrest cushions or a padded backrest cover.



Features and Options



Tension adjustable backrest

Adjustable in height and as wide as the seat.

Table tray

Ideal for use at meal and activity times.

Shorter chair

The chassis is 30 mm shorter, the seat depth is 50 mm shorter and the backrest 100 mm lower than the standard Rea Azalea. Shown here in the Electric Green frame colour.

A wide range of accessories can be found in our Accessories brochure, order using code 1576967.

Technical data


For more comprehensive pre-sales information about this product, including the product's user manual, please visit www.invacare.co.uk

Rea Azalea Minor	Seat width	Seat depth	Seat height	Backrest height	Armrest height	Length of legrests	Seat angle	Backrest angle
	340 - 440 mm	380 - 450 mm	400/450 mm	550 - 650 mm	240 - 360 mm	330 - 500 mm	+1° to 25°	0° - 30°
Rea Azalea Minor	SW + 220 mm (SW + 230 mm)	900 - 1250 mm	900 - 970 mm	32 kg	Max 75 kg	14,5 kg	ISO 7176-19	7°
	Driving unit width	Total height	Total length inc. Footrests	Total weight	Max. User weight	Transport weight	Crash tested	Maximum slope for use of parking brake

Frame colour



Upholstery



Invacare Limited
Pencoed Technology Park - Pencoed
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail: ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

