



INCREASED  
CLINICAL FUNCTIONALITY

# Invacare<sup>®</sup> Rea<sup>®</sup> Flex 3



Confidence in posture solutions

Trusted Solutions, Passionate People



Yes, you can.<sup>®</sup>

# FLEX 3

## Discover the Flex 3 system

The Flex 3 system has been specially designed to **fit a variety of people**, allowing a seat to be **tailored** to different body shapes and disabilities. Flex 3 aims to **improve posture and stability** which in turn increases activity levels and quality of life.





## Compatibility between wheelchairs

With the same attachment on all Invacare passive wheelchairs\*, the Flex 3 system provides the **flexibility** to combine backrest cushions with backrest types across its range. This allows for **backward compatibility** and **enhances refurbishment**.



\* Excluding Clematis® and Azalea® Minor



# Find the right wheel

## Rea Azalea® Range

The **Rea Azalea** range is the natural choice for the moderately active client as well as for the more passive and offers a **25° tilt** with a **30° low shear recline** function. The range offers a solution for a wide variety of requirements.

### **Azalea**

### **Azalea Base**

### **Azalea Minor**

### **Azalea Tall**

### **Azalea Max**

Adapts to a wide range of seating and positioning needs.

Accommodates alternative seating solutions.

Smaller dimensions for teenagers and petite adults.

Higher seat height and longer seat depths for taller people.

Meets the specific needs of the larger person.



# chair

## Rea Dahlia® Range

The **Rea Dahlia** is manoeuvrable, slim and adjustable with impressive tilt options and a **30° low shear recline** function. This versatile wheelchair is perfect for those seeking mobility with comfort and pressure relief.

**Dahlia 45** - Impressive **45° tilt** delivers excellent pressure relieving properties and supports the most vulnerable areas. To obtain pressure relief using only the tilt function, a wheelchair must go to the full 45° angle.

**Dahlia 30** - **30° tilt** version with lower floor to seat height for people who can readily propel themselves. With a self-tilt function available, clients can actively change posture themselves when required.



➡ **Now select the backrest system**

# Easy backrest guide



**Cover Lateral**  
The Cover Lateral follows the contoured shape set by the tension adjustable straps to give the client a good level of support.



**Laguna**  
The Laguna backrest cushion offers moderate lateral support and gives the client stability and support.



**Mistral 2**  
The Mistral 2 backrest cushion with its deep, reinforced wings and a lower, shaped upper section offers the client a distinct waist support and assists shoulder movement.



**Backrest Tension Adjustable**  
The backrest can be tailored to an individual shape by adjusting the straps. It is angle, width and height adjustable to provide a high level of support and comfort.



**Backrest Plate**  
A width and height adjustable, firm and stable support plate.





**Passad 2**  
 With its high shoulder support, the Passad 2 offers support for client with a weak posture. The high wings are reinforced.

**Shoulder High**  
 The Shoulder High backrest cushion is the same as Passad 2 but is 100 mm higher, and has been designed specially for taller clients.



## Proposed backrest options

Plate



**Short**

Plate



**Tall**

Either plate or tension adj



**Average**

Plate



**Wide**

Tension adj



**Wide above**

Tension adj



**Wide below**

Tension adj



**Larger in front**

Tension adj



**Larger behind**

Tension adj



**Kyphosis back**

Tension adj



**S-shaped back**

Take a look at our Flex 3  
video on our Invacare UK  
YouTube channel called  
**Invacare Flex 3 Seating System.**



**If further postural or seating support is needed, we  
have a wide range of accessories available. Take a look  
at our Rea accessory brochure for more information.**



Invacare Limited  
Pencoed Technology Park - Pencoed - Bridgend  
CF35 5AQ - United Kingdom  
Tel: +44 1 656 776 222  
Fax: +44 1656 776 220  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
Sales Order E-mail:  
[ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS - County Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

© 2016 Invacare International Sàrl  
All rights reserved.

All information quoted is believed  
to be correct at time of print.  
Invacare reserves the right to alter product  
specifications without prior consultation.  
Flex 3 - UK - 11/2016 - 1553824

Trusted Solutions,  
Passionate People



**Yes, you can.**