



Invacare[®] | **Comet[®] Ultra**

Built to support with performance
and reliability



Ideal for those with larger body shapes, the new **Comet^{Ultra}** offers a secure, supportive and powerful driving experience without compromising on style or functionality. Packed with all of the safety features and driving characteristics seen throughout the **Comet** family, the **Comet^{Ultra}** will get you to your destination quickly and efficiently, in complete comfort and control.

Features and options



Improved suspension

New multi adjustable suspension to offer superior comfort to larger individuals.



Reinforced seat column

Offers additional stability and safety when seated

Reinforced height and width adjustable armrest

Can be easily adjusted to the individual and offers additional support when getting on and off the scooter

Power, reliability and customisation

















Due to a carefully thought out design process and a robust construction, the safety and stability of the **Comet^{Ultra}** is never a concern; even for persons weighing up to 220kg. Developed with the highest quality components and safety features, the **Comet^{Ultra}** makes a trustworthy companion no matter where you go. Flexibility and adaptability means you can customise your scooter to suit your unique style and individual preference.

Complete comfort and support

The special contoured Ultra Seat has been developed in conjunction with therapists to ensure an optimum seating position and increased comfort. In addition, a pneumatically adjustable lumbar support allows for further comfort when seated if and when required. The reinforced seat base and adjustable armrests offer excellent support and allows the use of the armrest as additional support when getting on and off the scooter.

Technical data

For more comprehensive information about this product, including the product's user manual, please visit www.invacare.co.uk

	Backrest width	Seat depth	Seat to floor height	Backrest height	Driving unit width	Total length	Empty weight	The heaviest part
Comet^{Ultra}	 660 mm	 535 mm	 Min. 440 mm Max. 510 mm	 515 mm	 665 - 685 mm	 1450 mm	 148 kg	 approx. 63 kg
	Max. user weight	Battery capacity	Turning diameter	Max. climbable obstacle height	Ground clearance	Max. safe slope	Max. driving distance	Max. driving speed
Comet^{Ultra}	 220 kg	 2 x 12 V/75 Ah (C20) AGM	 2750 mm	 100 mm	 100 mm	 10 °	 58 km (36 miles)	 10 km/h (6mph)

Motor capacity



Comet^{Ultra} 10km/h: 650 W / 1800 W¹

1. Peak performance

Colour of shrouds

Classic Colour Line

Marble White Indian Blue Jade Green Diamond Silver Jasper Red



Premium Colour Line

Carmine Red Iconic Chrome Electric Blue Luxor Beige Crystal Copper



Upholstery

Black Vinyl



Invacare Limited
Pencoed Technology Park - Pencoed
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS -
County Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

Trusted Solutions,
Passionate People



© 2017 Invacare International Sàrl. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Comet Ultra - UK - 03/2017 - 1627900