



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

APRIL 2020
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Kuschall Champion Aluminium

Standard Features

- Max user weight 120kg/19 Stone
- Active, lightweight, folding aluminium wheelchair
- Multiple adjustments, front & rear seat heights and centre of gravity
- Aluminium front frame with 75° or 85° front angle options
- 24" Self propel Q/R wheels with grey aluminium handrims
- Schwalbe black & grey pneumatic tyres
- Choice of 7 seat widths 360mm to 480mm (14" - 19")
- Choice of 8 seat depths 340mm to 480mm (13" - 19")
- Choice of 13 frame colours
- Range of front seat heights from 450mm to 540mm (17" - 21")
- Range of rear seat heights from 390mm to 490mm (15" - 19")
- Angle adjustable, soft locking, folding backrest with nylon light upholstery
- Black comfort front castor
- Standard brake lever bent
- Tool kit



Standard Base price: £2,220

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
DDC0031	KÜSCHALL CHAMPION	£2,220
DDC0000	MYKÜSCHALL SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE	POA

FRAME*				
DDC0203		FRAME ANGLE 75° ALUMINIUM	(NC 0138)	◆
DDC0212		FRAME ANGLE 85° ALUMINIUM	(NC 0138)	○
DDC0201		FRAME ANGLE 75° ALUMINIUM - TAPERED FRONT FRAME	(NC 0126, 0128, 1722, 1714)	○
DDC0243		FRAME ANGLE 85° ALUMINIUM - TAPERED FRONT FRAME	(NC 0126, 1722, 1714)	○
SEAT WIDTH*				
DDC0126		360 MM	(NC 1722)	○
DDC0128		380 MM	(NC 1722)	○
DDC0130		400 MM		○
DDC0132		420 MM		○
DDC0134		440 MM		○
DDC0136		460 MM		○
DDC0138		480 MM		○
SEAT DEPTH*				
DDC0410		340 MM	(NC 0615, 0662, 0663)	○
DDC0412		360 MM	(NC 0662, 0663)	○
DDC0415		380 MM	(NC 0662, 0663)	○
DDC0417		400 MM		○
DDC0419		420 MM		○
DDC0421		440 MM		○
DDC0424		460 MM		○
DDC0427		480 MM		○
BACK HEIGHT*				
DDC0585		BACK HEIGHT FOR FIXED BACK FINAL HEIGHT DEPENDS ON THE HEIGHT OF THE POSITIONING BACKREST, SELECT IF MATRX BACKREST IS REQUIRED.	(NC 0700, 0726, 0735)	○
DDC0550		300 MM	(NC 0700, 0769, 0750, 1491)	○
DDC0503		315 MM	(NC 0700, 0769, 0750, 1491)	○
DDC0506		330 MM	(NC 0769, 1491)	○
DDC0509		345 MM	(NC 0769, 1491)	○
DDC0512		360 MM	(NC 0769, 1491)	○
DDC0514		375 MM	(NC 0769, 1491)	○
DDC0516		390 MM	(NC 0769, 1491)	○
DDC0518		405 MM	(NC 0769)	○
DDC0520		420 MM	(NC 0769)	○
DDC0522		435 MM	(NC 0769)	○
DDC0524		450 MM	(NC 0769)	○
DDC0526		465 MM	(NC 0769)	○

FRONT SEAT HEIGHT*				
DDC0847		450 MM	(NC 1684, 1685, 1686, 1687, 1688, 1689, 1690, 1691)	<input type="radio"/>
DDC0848		460 MM	(NC 1685, 1686, 1687, 1688, 1689, 1690, 1691)	<input type="radio"/>
DDC0849		470 MM	(NC 1686, 1687, 1688, 1689, 1690, 1691)	<input type="radio"/>
DDC0850		480 MM	(NC 1687, 1688, 1689, 1690, 1691)	<input type="radio"/>
DDC0851		490 MM	(NC 1688, 1689, 1690, 1691)	<input type="radio"/>
DDC0852		500 MM	(NC 1689, 1690, 1691)	<input type="radio"/>
DDC0853		510 MM	(NC 1690, 1691)	<input type="radio"/>
DDC0854		520 MM	(NC 1691)	<input type="radio"/>
DDC0855		530 MM	(NC 0800)	<input type="radio"/>
DDC0856		540 MM POSSIBLE WITH 5" FRONT WHEELS ONLY	(NC 0800, 0801, 0860, 0809, 0861, 0803, 0808, 0825, 0807, 0831, 0796, 0797)	<input type="radio"/>
REAR SEAT HEIGHT*				
DDC0944		390 MM POSSIBLE WITH 22" & 24" REAR WHEELS	(NC 0925, 0903)	<input type="radio"/>
DDC0945		400 MM POSSIBLE WITH 22" & 24" REAR WHEELS	(NC 0925, 0903)	<input type="radio"/>
DDC0946		410 MM POSSIBLE WITH 22", 24" & 25" REAR WHEELS	(NC 0903)	<input type="radio"/>
DDC0947		420 MM POSSIBLE WITH 22", 24" & 25" REAR WHEELS	(NC 0903)	<input type="radio"/>
DDC0948		430 MM		<input type="radio"/>
DDC0949		440 MM		<input type="radio"/>
DDC0950		450 MM POSSIBLE WITH 24", 25" & 26" REAR WHEELS	(NC 0901)	<input type="radio"/>
DDC0951		460 MM POSSIBLE WITH 24", 25" & 26" REAR WHEELS	(NC 0901)	<input type="radio"/>
DDC0952		470 MM POSSIBLE WITH 25" & 26" REAR WHEELS	(NC 0901, 0902)	<input type="radio"/>
DDC0953		480 MM POSSIBLE WITH 25" & 26" REAR WHEELS	(NC 0901, 0902)	<input type="radio"/>
DDC0954		490 MM POSSIBLE WITH 26" REAR WHEELS ONLY	(NC 0901, 0902, 0925)	<input type="radio"/>
LOWER LEG LENGTH*				
DDC1673		320 MM		<input type="radio"/>
DDC1674		330 MM		<input type="radio"/>
DDC1675		340 MM		<input type="radio"/>
DDC1676		350 MM		<input type="radio"/>
DDC1677		360 MM		<input type="radio"/>
DDC1678		370 MM		<input type="radio"/>
DDC1679		380 MM		<input type="radio"/>
DDC1680		390 MM		<input type="radio"/>
DDC1681		400 MM		<input type="radio"/>
DDC1682		410 MM		<input type="radio"/>
DDC1683		420 MM		<input type="radio"/>
DDC1684		430 MM		<input type="radio"/>
DDC1685		440 MM		<input type="radio"/>
DDC1686		450 MM		<input type="radio"/>
DDC1687		460 MM		<input type="radio"/>
DDC1688		470 MM		<input type="radio"/>
DDC1689		480 MM		<input type="radio"/>
DDC1690		490 MM		<input type="radio"/>
DDC1691		500 MM		<input type="radio"/>

BACKREST*			
DDC0734		ANGLE ADJUSTABLE, FORWARD FOLDING WITH SOFT LOCK. BACKREST REMAINS LOCKED FOR TRANSPORTATION	◆
BACKREST UPHOLSTERY*			
DDC0700		BLACK STANDARD LIGHT UPHOLSTERY, WITHOUT TENSION ADJUSTMENT (NC 0703, 1491)	◆
DDC0726		TENSION ADJUSTABLE BACKREST, WITH PREMIUM BLACK NYLON AND REFLECTIVE TRIM UPHOLSTERY COVER	£113
DDC0735		TENSION ADJUSTABLE BACKREST, WITH BLACK AERO/SPACER UPHOLSTERY COVER	£185
DDC0769		WITHOUT BACKREST UPHOLSTERY SELECT IF ORDERING MATRX BACKREST, OR IF NO UPHOLSTERY IS REQUIRED (NC 0751, 0702, 0703, 0704, 1491)	○
PUSH HANDLES*			
DDC0750		WITHOUT PUSH HANDLES	◆
DDC0751		PUSH HANDLES FIXED SHORT/MINI HAND GRIP	£110
DDC0702		PUSH HANDLES FIXED LONG HAND GRIP	○
DDC0703		PUSH HANDLES HEIGHT ADJUSTABLE, REAR SET (NC 1411)	£210
DDC0704		PUSH HANDLES FOLDABLE (NC 1491)	£164
DDC0733		PUSH HANDLES HEIGHT ADJUSTABLE, INTEGRATED (NC 1491)	£181
BACKREST ANGLE*			
DDC0713		76°	○
DDC0714		80.5°	○
DDC0715		85°	◆
DDC0716		89.5°	○
DDC0717		94°	○
BACKREST STABILIZATION BAR			
DDC1491		BACKREST STABILIZATION BAR, SWINGAWAY AND REMOVABLE	£136
SEAT UPHOLSTERY*			
DDC0610		STANDARD SEAT UPHOLSTERY	◆
DDC0615		LIGHT SEAT UPHOLSTERY, MESH MATERIAL TO REDUCE WEIGHT	£57
SEAT CUSHION			
DDC0600		STANDARD SEAT CUSHION, 60 MM DOUBLE FOAM WITH VISCO TOP LAYER	£55
DDC0662		MATRX LIBRA CUSHION, BLACK SMOOTH COVER - VERY HIGH RISK	£323
DDC0663		MATRX VI CUSHION, BLACK SMOOTH COVER - HIGH RISK	£299
SEAT BAG			
DDC0606		SEAT BAG, CLICK IN AND REMOVABLE, INCLUDES THE FIXATION CLIPS (NC 1714)	£61
FOOTREST*			
DDC1500		RIGID, PLASTIC FOOTPLATE, HEIGHT AND ANGLE ADJUSTABLE, SET TO THE FRONT. NOTE: THIS DOES NOT FLIP UP FOR TRANSFER, BUT IT FOLDS EASY FOR TRANSPORT AND STORAGE	○
DDC1501		RIGID, PLASTIC FOOTPLATE, HEIGHT AND ANGLE ADJUSTABLE, SET TO THE REAR NOTE: THIS DOES NOT FLIP UP FOR TRANSFER, BUT IT FOLDS EASY FOR TRANSPORT AND STORAGE	◆
CALF STRAP*			
DDC1533		CALF STRAP MOUNTED ON SUPPORTER	◆
DDC1535		CALF STRAP MOUNTED ON FRAME	○

CLOTHESGUARDS/MUDGUARDS*			
DDC1401		CARBON CLOTHES GUARDS FLIP BACK	◆
DDC1405		CARBON MUDGUARDS FIXED (NC 1416)	£205
DDC1408		CARBON MUDGUARDS REMOVABLE (NC 1416)	£255
DDC1111		NO CLOTHESGUARDS SELECT WITH OPTIONS 1416 ARMREST	○
ARMREST			
DDC1411		TUBULAR ARMREST. HEIGHT ADJUSTABLE, REMOVABLE, NON-LOCKING, SWING-AWAY AND COMFORT PADDED, NO CLOTHESGUARD INCLUDED (NC 1111)	£199
DDC1416		ARMREST, HEIGHT ADJUSTABLE, SHORT ARMPAD (T-ARM), NON-LOCKING, NO CLOTHESGUARD INCLUDED	£88
FRAME COLOUR*			
DDC0300		GLOSSY-BLACK	○
DDC0302		ROYAL BLUE	○
DDC0303		ROYAL RED	○
DDC0314		BLACK MATT	○
DDC0315		SILVER GREY	○
DDC0366		DYNAMITE ORANGE	○
DDC0365		MERCURY BLUE	○
DDC0368		NITRO PINK	○
DDC0367		ATOMIC VIOLET	○
DDC0369		MINT GREEN	○
DDC0356		IVORY WHITE	○
DDC0358		SMOKEY GREY	○
DDC0359		FLASHY YELLOW	○
COLOUR POWER SET			
DDC0338		CHROME STAREC REAR WHEEL RIM ONLY AVAILABLE IN BLACK WHEN CPS CHROME IS SELECTED	£185
DDC0330		TECHNO-ANTHRACITE	£185
DDC0329		COMET RED	£185
DDC0339		CHAMPAGNE	£185
DDC0334		B-BLUE	£185
DDC0336		ARCTIC BLUE	£185

FRONT CASTORS*			
DDC0800		75MM X 25MM (3" X 1") SUPERSPORT CASTOR	○
DDC0801		100MM X 30MM (4" X 1.2") EVERYDAY CASTOR WITH BLACK TYRE	○
DDC0860		100MM X 30MM (4" X 1.2") COMFORT CASTOR WITH BLACK TYRE	○
DDC0809		100MM X 25MM (4" X 1") STAREC CASTOR WITH BLACK HUB, BLACK TYRE	£69
DDC0861		100MM X 25MM (4" X 1") STAREC CASTOR WITH CHROME BLACK HUB, BLACK TYRE	£69
DDC0803		100MM X 25MM (4" X 1") STAREC CASTOR WITH TECHNO-ANTHRACITE BLACK HUB, BLACK TYRE	£69
DDC0808		100MM X 25MM (4" X 1") STAREC CASTOR WITH COMET RED BLACK HUB, BLACK TYRE	£69
DDC0825		100MM X 25MM (4" X 1") STAREC CASTOR WITH CHAMPAGNE BLACK HUB, BLACK TYRE	£69
DDC0807		100MM X 25MM (4" X 1") STAREC CASTOR WITH B-BLUE BLACK HUB, BLACK TYRE	£69
DDC0831		100MM X 25MM (4" X 1") STAREC CASTOR WITH ARCTIC BLUE BLACK HUB, BLACK TYRE	£69
DDC0796		100MM X 36MM (4" X 1.4") ALUMINIUM SOFTROLL CASTOR WITH BLACK HUB, BLACK TYRE	£77
DDC0797		100MM X 36MM (4" X 1.4") ALUMINIUM SOFTROLL CASTOR WITH SILVER HUB, BLACK TYRE	£77
DDC0810		125MM X 30MM (5" X 1.2") EVERYDAY CASTOR WITH BLACK TYRE	○
DDC0863		125MM X 30MM (5" X 1.2") COMFORT CASTOR WITH BLACK TYRE	◆
DDC0819		125MM X 25MM (5" X 1") STAREC CASTOR BLACK HUB, BLACK TYRE	£69
DDC0862		125MM X 25MM (5" X 1") STAREC CASTOR CHROME HUB, BLACK TYRE	£69
DDC0813		125MM X 25MM (5" X 1") STAREC CASTOR TECHNO-ANTHRACITE HUB, BLACK TYRE	£69
DDC0818		125MM X 25MM (5" X 1") STAREC CASTOR COMET RED HUB, BLACK TYRE	£69
DDC0827		125MM X 25MM (5" X 1") STAREC CASTOR CHAMPAGNE HUB, BLACK TYRE	£69
DDC0817		125MM X 25MM (5" X 1") STAREC CASTOR B-BLUE HUB, BLACK TYRE	£69
DDC0832		125MM X 25MM (5" X 1") STAREC CASTOR ARCTIC BLUE HUB, BLACK TYRE	£69
DDC0798		125MM X 36MM (5" X 1.4") ALUMINIUM SOFTROLL CASTOR WITH BLACK HUB, BLACK TYRE	£77
DDC0799		125MM X 36MM (5" X 1.4") ALUMINIUM SOFTROLL CASTOR WITH SILVER HUB, BLACK TYRE	£77
REAR WHEEL POSITION*			
DDC0263		1ST FROM FRONT MOST ACTIVE POSITION	○
DDC0264		2ND FROM FRONT	○
DDC0265		3RD FROM FRONT	○
DDC0266		4TH FROM FRONT	○
DDC0267		5TH FROM FRONT	◆
DDC0268		6TH FROM FRONT	○
DDC0269		7TH FROM FRONT	○
DDC0259		8TH FROM FRONT	○
DDC0260		9TH FROM FRONT LEAST ACTIVE POSITION	○
QUAD QUICK RELEASE AXLE			
DDC0909		QUAD QUICK RELEASE/TETRA CLIP	(NC 0906, 0912) £71

REAR WHEEL SIZE*				
DDC0901		22"	(NC 0908, 0918, 0906, 0912, 1010, 1004, 1011, 1104, 1106, 1121, 1122, 1202, 1208)	○
DDC0902		24"	(NC 0918, 0912)	◆
DDC0925		25"	(NC 0908, 0918, 0906, 1002, 1004, 1102, 1104, 1106, 1121, 1122)	○
DDC0903		26"	(NC 0908, 0918, 0906, 0912, 1002, 1004, 1011, 1106, 1121, 1122, 1202, 1208)	○
REAR WHEEL CAMBER*				
DDC0282		0° CAMBER		◆
DDC0279		2° CAMBER		○
DDC0283		4° CAMBER	(NC 1302)	○
REAR WHEEL TYPE*				
DDC0915		22", 24", 26" KÜSCHALL STAREC, CROSS SPOKED STAREC HUB COLOURLESS AND RIM BLACK - HUB AND RIM IN CPS WHEN SELECTED		◆
DDC0908		24" HIGH PERFORMANCE WHEEL STAREC HUB AND RIM BLACK - OR IN CPS WHEN SELECTED. NOTE: RIM SUPPLIED IN BLACK WITH CHROME CPS	(NC 1004)	£210
DDC0918		25" HIGH PERFORMANCE WHEEL STAREC HUB AND RIM BLACK - HUB IN CPS WHEN SELECTED, BLACK RIM	(NC 1004, 1200, 1202, 1208)	£302
DDC0906		24" SPINERGY LXX HUB, RIM AND SPOKES BLACK - CPS NOT APPLICABLE IF CHOSEN	(NC 1004, 1202, 1208)	£505
DDC0912		25" SPINERGY LXX HUB, RIM, SPOKES BLACK - CPS NOT APPLICABLE IF CHOSEN	(NC 1004, 1200, 1202, 1208)	£505
SWO0907		24" OR 25" SPINERGY SPOX WHEEL (CHROME HUB, 18 SPOKES IN BLACK, NO CPS) COLOURED SPOKES ARE AVAILABLE VIA MYKUSCHALL, PLEASE INDICATE REQUIREMENTS IN THE SPACE PROVIDED AT THE END OF THE FORM	(NC 1004, 1200, 1202, 1208)	£402
SWO0909		24" OR 25" SPINERGY LX WHEEL (CHROME HUB, 12 SPOKES IN BLACK, NO CPS) COLOURED SPOKES ARE AVAILABLE VIA MYKUSCHALL, PLEASE INDICATE REQUIREMENTS IN THE SPACE PROVIDED AT THE END OF THE FORM	(NC 1004, 1200, 1202, 1208)	£608
REAR WHEEL TYRE*				
DDC1014		22", 24", 25", 26" SCHWALBE RIGHTRUN BLACK/GREY PNEUMATIC TYRE		○
DDC1002		22", 24" PROFILE GREY PNEUMATIC TYRE 7 BAR		◆
DDC1010		24", 25", 26" SCHWALBE MARATHON PLUS SMART GUARD (PUNCTURE PROTECTED TYRE WITH 5MM LAYER OF FLEXIBLE RUBBER)		£82
DDC1004		24" SOLID PROFILE GREY TYRE		£55
DDC1011		24", 25" SCHWALBE ONE V-GUARD PUNCTURE RESISTANT TYRE (CUT-RESISTANT HIGH-TECH FIBRE)		£80

HANDRIM*			
DDC1102		22", 24", 26" SUPERGRIPP, MOUNTED AT 3 CM DISTANCE	£185
DDC1104		24", 26" STEEL, MOUNTED AT 3 CM DISTANCE	£52
DDC1106		24" TITANIUM, MOUNTED AT 3 CM DISTANCE	£159
DDC1121		24" CARBOLIFE GEKKO MOUNTED AT 3 CM DISTANCE	£345
DDC1113		22", 24", 25", 26" LIGHT ALUMINIUM GREY, MOUNTED AT 3 CM DISTANCE	◆
DDC1117		22", 24", 25", 26" LIGHT ALUMINIUM BLACK, MOUNTED AT 3 CM DISTANCE	£34
DDC1114		22", 24", 25", 26" LIGHT ALUMINIUM GREY, MOUNTED AT 4 CM DISTANCE	○
DDC1122		24" CARBOLIFE CURVE L, MOUNTED AT 3 CM DISTANCE	£308
SPOKEGUARDS			
DDC1200		22", 24", 26" TRANSPARENT	£104
DDC1202		24" KÜSCH	£104
DDC1208		24" TRANSPARENT AND WHITE "KÜSCHALL"	£104
FIXATION CLIP			
DDC1722		FIXATION STRAP (HOOK AND LOOP STRAP UNDER THE SEAT HOLDS THE FRAME TOGETHER WHEN FOLDED FOR EASY TRANSPORTATION, DESIGNED FOR CLIENTS WITH LIMITED HAND FUNCTION)	○
DDC1714		FIXATION CLIP (PLASTIC CLIP ON EITHER SIDE OF THE SEAT HOLDS THE FRAME TOGETHER WHEN FOLDED FOR EASY TRANSPORTATION)	○
UNFOLDING HELP*			
DDC1719		UNFOLDING HELP (CORD TO HELP OPEN THE FOLDING MECHANISM FROM THE FRONT OF THE CHAIR)	◆
BRAKES*			
DDC1300		STANDARD PUSH TO LOCK BRAKES, LEVER BENT	◆
DDC1301		STANDARD PUSH TO LOCK BRAKES, LEVER STRAIGHT	○
DDC1306		PULL BRAKES	£64
DDC1302		PERFORMANCE BRAKES, PUSH TO LOCK. BRAKE BARREL FLIPS OUT OF THE WAY OF WHEEL (NC 1310)	£64
DDC1303		ACTIVE/SCISSOR BRAKES (NC 1310)	£64
BRAKE LEVER EXTENSION			
DDC1310		BRAKE LEVER EXTENSION	£41
SAFETY OPTIONS			
DDC1702		ANTI-TIPPER RIGHT	£101
DDC1703		ANTI-TIPPER LEFT	£101
DDC1704		ANTI-TIPPER BOTH SIDES	£202
DDC1730		STANDARD PELVIC BELT	£61

OTHER OPTIONS			
DDC1739		FRONT FRAME PROTECTORS; 2 PIECE. REMOVABLE WITH HOOK AND LOOP	£31
DDC1708		TIPPER-AID, RIGHT	£38
DDC1706		CANE HOLDER, MOUNTED RIGHT	£62
DDC1715		AIR PUMP	£11
DDC1712		KÜSCHALL WATCH, RIGHT SIDE	£88
DDC1713		KÜSCHALL WATCH, LEFT SIDE	£88
DDC1717		TOOL KIT	◆
DDC1728		WHEEL BAG (FITS ALL WHEEL SIZES)	£45
DDC1734		24" MOUNTAIN BIKE WHEELS WITH BAG. FITTED WITH SPECIFIC BLACK TYRES AND GREY ALUMINIUM HANDRIMS ; QUICK RELEASE AXLES INCLUDED	£265
DDC1740		TRAVEL BAG FOR CHAMPION SOFT BAG, WHEELS CAN REMAIN ON THE WHEELCHAIR	£210

TOTAL PRICE	
--------------------	--

Non Conforming Options:

DDC0126	360 MM
DDC0128	380 MM
DDC0138	480 MM
DDC1684	430 MM
DDC1685	440 MM
DDC1686	450 MM
DDC1687	460 MM
DDC1688	470 MM
DDC1689	480 MM
DDC1690	490 MM
DDC1691	500 MM
DDC0700	BLACK STANDARD LIGHT UPHOLSTERY, WITHOUT TENSION ADJUSTMENT
DDC0726	TENSION ADJUSTABLE BACKREST, WITH PREMIUM BLACK NYLON AND REFLECTIVE TRIM UPHOLSTERY COVER
DDC0735	TENSION ADJUSTABLE BACKREST, WITH BLACK AERO/SPACER UPHOLSTERY COVER
DDC0769	WITHOUT BACKREST UPHOLSTERY SELECT IF ORDERING MATRX BACKREST, OR IF NO UPHOLSTERY IS REQUIRED
DDC0750	WITHOUT PUSH HANDLES
DDC0751	PUSH HANDLES FIXED SHORT/MINI HAND GRIP
DDC0702	PUSH HANDLES FIXED LONG HAND GRIP
DDC0703	PUSH HANDLES HEIGHT ADJUSTABLE, REAR SET
DDC0704	PUSH HANDLES FOLDABLE
DDC1491	BACKREST STABILIZATION BAR, SWINGAWAY AND REMOVABLE
DDC0615	LIGHT SEAT UPHOLSTERY, MESH MATERIAL TO REDUCE WEIGHT
DDC0662	MATRX LIBRA CUSHION, BLACK SMOOTH COVER - VERY HIGH RISK
DDC0663	MATRX VI CUSHION, BLACK SMOOTH COVER - HIGH RISK
DDC1111	NO CLOTHESGUARDS SELECT WITH OPTIONS 1416 ARMREST
DDC1411	TUBULAR ARMREST. HEIGHT ADJUSTABLE, REMOVABLE, NON-LOCKING, SWING-AWAY AND COMFORT PADDED, NO CLOTHESGUARD INCLUDED
DDC1416	ARMREST, HEIGHT ADJUSTABLE, SHORT ARMPAD (T-ARM), NON-LOCKING, NO CLOTHESGUARD INCLUDED
DDC0800	75MM X 25MM (3" X 1") SUPERSPORT CASTOR
DDC0801	100MM X 30MM (4" X 1.2") EVERYDAY CASTOR WITH BLACK TYRE
DDC0860	100MM X 30MM (4" X 1.2") COMFORT CASTOR WITH BLACK TYRE
DDC0809	100MM X 25MM (4" X 1") STAREC CASTOR WITH BLACK HUB, BLACK TYRE
DDC0861	100MM X 25MM (4" X 1") STAREC CASTOR WITH CHROME BLACK HUB, BLACK TYRE
DDC0803	100MM X 25MM (4" X 1") STAREC CASTOR WITH TECHNO-ANTHRACITE BLACK HUB, BLACK TYRE
DDC0808	100MM X 25MM (4" X 1") STAREC CASTOR WITH COMET RED BLACK HUB, BLACK TYRE
DDC0825	100MM X 25MM (4" X 1") STAREC CASTOR WITH CHAMPAGNE BLACK HUB, BLACK TYRE
DDC0807	100MM X 25MM (4" X 1") STAREC CASTOR WITH B-BLUE BLACK HUB, BLACK TYRE
DDC0831	100MM X 25MM (4" X 1") STAREC CASTOR WITH ARCTIC BLUE BLACK HUB, BLACK TYRE
DDC0796	100MM X 36MM (4" X 1.4") ALUMINIUM SOFTROLL CASTOR WITH BLACK HUB, BLACK TYRE
DDC0797	100MM X 36MM (4" X 1.4") ALUMINIUM SOFTROLL CASTOR WITH SILVER HUB, BLACK TYRE
DDC0901	22"
DDC0902	24"
DDC0925	25"

DDC0903	26"
DDC0908	24" HIGH PERFORMANCE WHEEL STAREC HUB AND RIM BLACK - OR IN CPS WHEN SELECTED. NOTE: RIM SUPPLIED IN BLACK WITH CHROME CPS
DDC0918	25" HIGH PERFORMANCE WHEEL STAREC HUB AND RIM BLACK - HUB IN CPS WHEN SELECTED, BLACK RIM
DDC0906	24" SPINERGY LXX HUB, RIM AND SPOKES BLACK - CPS NOT APPLICABLE IF CHOSEN
DDC0912	25" SPINERGY LXX HUB, RIM, SPOKES BLACK - CPS NOT APPLICABLE IF CHOSEN
DDC1002	22", 24" PROFILE GREY PNEUMATIC TYRE 7 BAR
DDC1010	24", 25", 26" SCHWALBE MARATHON PLUS SMART GUARD (PUNCTURE PROTECTED TYRE WITH 5MM LAYER OF FLEXIBLE RUBBER)
DDC1004	24" SOLID PROFILE GREY TYRE
DDC1011	24", 25" SCHWALBE ONE V-GUARD PUNCTURE RESISTANT TYRE (CUT-RESISTANT HIGH-TECH FIBRE)
DDC1102	22", 24", 26" SUPERGRIPP, MOUNTED AT 3 CM DISTANCE
DDC1104	24", 26" STEEL, MOUNTED AT 3 CM DISTANCE
DDC1106	24" TITANIUM, MOUNTED AT 3 CM DISTANCE
DDC1121	24" CARBOLIFE GEKKO MOUNTED AT 3 CM DISTANCE
DDC1122	24" CARBOLIFE CURVE L, MOUNTED AT 3 CM DISTANCE
DDC1200	22", 24", 26" TRANSPARENT
DDC1202	24" KÜSCH
DDC1208	24" TRANSPARENT AND WHITE "KÜSCHALL"
DDC1722	FIXATION STRAP (HOOK AND LOOP STRAP UNDER THE SEAT HOLDS THE FRAME TOGETHER WHEN FOLDED FOR EASY TRANSPORTATION, DESIGNED FOR CLIENTS WITH LIMITED HAND FUNCTION)
DDC1714	FIXATION CLIP (PLASTIC CLIP ON EITHER SIDE OF THE SEAT HOLDS THE FRAME TOGETHER WHEN FOLDED FOR EASY TRANSPORTATION)
DDC1302	PERFORMANCE BRAKES, PUSH TO LOCK. BRAKE BARREL FLIPS OUT OF THE WAY OF WHEEL
DDC1310	BRAKE LEVER EXTENSION



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):