


Tel: 01656 776222 Fax: 01656 776220  
Email: ordersuk@invacare.com  
www.invacare.co.uk

<b>Customer name:</b>	<b>Account No :</b>
<b>Contact no:</b>	<b>Order No :</b>
<b>Delivery address:</b>	<b>Order Qty :</b>
	<b>Order Date :</b>

**Matrx Libra**  
Retail order Form  
Maximum User Weight 21 stone (136kg) / HD 36 stone (227kg)


		WIDTH								
cm		36	38	41	43	46	48	51	53	56
inch		14	15	16	17	18	19	20	21	22
DEPTH	30									
	36									
	41									
	43									
	46									
	51									
56										



**Ordering the correct size**

Please note the complete product dimensions. The example shown is of a 16"x16" cushion part number LC1616. Please allow an additional 1 1/2" to the required seat depth for the rear cushion radius.

The rear cushion radius is a curved extension designed to maximise surface area and minimise the gap between seat and backrest; while stabilising fluid in tilted positions, with coccyx cut out for added relief.



- Features:**
- Very High Risk
  - 3" height
  - 3 Year Warranty
  - Anatomically contoured foam base
  - Innovative dual layer fluid sac
  - Choice of 3 covers
  - Positioning accessories available



Code	X	Description	Price
<b>Libra Cushion</b>			
LC _____ (Width x Depth)		Matrx Libra standard sizes with black smooth cover	£354
LCW _____ (Width x Depth)		Matrx Libra standard sizes with black wipedown cover	£354
LCS _____ (Width x Depth)		Matrx Libra standard sizes with black spacer cover	£354

## Heavy Duty sizes

LC2117	Heavy Duty Libra 21x17" (53x43cm)	Black Smooth cover	£486
LC2118	Heavy Duty Libra 21x18" (53x46cm)	Black Smooth cover	£486
LC2120	Heavy Duty Libra 21x20" (53x51cm)	Black Smooth cover	£486
LC2122	Heavy Duty Libra 21x22" (53x56cm)	Black Smooth cover	£486
LC2217	Heavy Duty Libra 22x17" (56x43cm)	Black Smooth cover	£486
LC2218	Heavy Duty Libra 22x18" (56x46cm)	Black Smooth cover	£486
LC2220	Heavy Duty Libra 22x20" (56x51cm)	Black Smooth cover	£486
LC2222	Heavy Duty Libra 22x22" (56x56cm)	Black Smooth cover	£486
LCW2117	Heavy Duty Libra 21x17" (53x43cm)	Black Wipedown Cover	£486
LCW2118	Heavy Duty Libra 21x18" (53x46cm)	Black Wipedown Cover	£486
LCW2120	Heavy Duty Libra 21x20" (53x51cm)	Black Wipedown Cover	£486
LCW2122	Heavy Duty Libra 21x22" (53x56cm)	Black Wipedown Cover	£486
LCW2217	Heavy Duty Libra 22x18" (56x43cm)	Black Wipedown Cover	£486
LCW2218	Heavy Duty Libra 22x18" (56x46cm)	Black Wipedown Cover	£486
LCW2220	Heavy Duty Libra 22x20" (56x51cm)	Black Wipedown Cover	£486
LCW2222	Heavy Duty Libra 22x22" (56x56cm)	Black Wipedown Cover	£486
LCS2117	Heavy Duty Libra 21x17" (53x43cm)	Spacer Fabric Cover	£486
LCS2118	Heavy Duty Libra 21x18" (53x46cm)	Spacer Fabric Cover	£486
LCS2120	Heavy Duty Libra 21x20" (53x51cm)	Spacer Fabric Cover	£486
LCS2122	Heavy Duty Libra 21x22" (53x56cm)	Spacer Fabric Cover	£486
LCS2217	Heavy Duty Libra 22x17" (56x43cm)	Spacer Fabric Cover	£486
LCS2218	Heavy Duty Libra 22x18" (56x46cm)	Spacer Fabric Cover	£486
LCS2220	Heavy Duty Libra 22x20" (56x51cm)	Spacer Fabric Cover	£486
LCS2222	Heavy Duty Libra 22x22" (56x56cm)	Spacer Fabric Cover	£486

## Accessories



LABD	<b>Matrix Abductor</b> For comfortable mid-line positioning of thighs and hips (Provides 1 1/8" (29mm) increased height abductor contour)		<b>£35</b>
1487157	<b>Anterior Seat Wedge</b> Increases slope, prevents sliding and reduces shear (1" Anterior wedge increases pre-ischial contour - Underside mounted)		<b>£26</b>
1487158	<b>Tapered Wedges (pair)</b> For obliquity and pelvic rotation when used in rear (3/4" Tapered wedge - Underside mounted)		<b>£18</b>
1487155	<b>Lateral Wedge</b> Reinforces lateral stability for enhanced midline positioning (1" Lateral wedge increases hip/thigh loading, - Underside mounted)		<b>£29</b>
1578357	<b>Libra Fluid Supplement</b> For targeted protection/obliquity accommodation (May be used to accommodate mild pelvic obliquity/additional pressure reduction)		<b>£48</b>
1589366	<b>Libra Accessory Kit</b> Comprises of all the above items in an assessment presentation box with instruction guide. Parts included: LABD, 1487157, 1487158, 1487155, 1578357		<b>£107</b>
	<b>Spare Plaitlon Inner Liner</b> Provides added protection to the foam cushion from infection control issues.		<b>£12</b>
	<b>Additional Cover (Widths 14" to 20")</b>	<input type="checkbox"/> Black Smooth <input type="checkbox"/> Black Spacer <input type="checkbox"/> Black Wipedown	<b>£91</b>
	<b>Additional Cover (Widths 21" and over )</b>	<input type="checkbox"/> Black Smooth <input type="checkbox"/> Black Spacer <input type="checkbox"/> Black Wipedown	<b>£107</b>
1494494	14" x 16" Rigidiser		<b>£37</b>
1487165	16" x 16" Rigidiser		<b>£37</b>
1494495	16" x 18" Rigidiser		<b>£37</b>
1494496	16" x 20" Rigidiser		<b>£72</b>
1494497	18" x 16" Rigidiser		<b>£37</b>
1494498	18" x 18" Rigidiser		<b>£37</b>
1494499	18" x 20" Rigidiser		<b>£72</b>
1494500	20" x 16" Rigidiser		<b>£37</b>
1494501	20" x 18" Rigidiser		<b>£37</b>
1494502	20" x 20" Rigidiser		<b>£72</b>
1494503	22" x 18" Rigidiser		<b>£72</b>
1494504	22" x 20" Rigidiser		<b>£72</b>
<b>TOTAL PRICE:</b>			

**Special Instructions - Extended lead times will apply**

**Name:**

**Email:**

**Phone:**



**Description of special request (please include relevant option code):**